

## Exhibit C

NOTEBOOK NO. 2369

ISSUED TO RANDY SAIKI

ON \_\_\_\_\_ 19\_\_\_\_

DEPARTMENT HUMAN GENETICS

RETURNED \_\_\_\_\_ 19 \_\_\_\_\_

— SCIENTIFIC NOTEBOOK CO. —

**5007 WEST DONNA DRIVE**

**STEVENSVILLE, MICHIGAN 49127**

From Page No. X

Rec'd 80 $\mu$ l Tag polymerase from David Gelfand. Tube is labeled "fraction VIII A", will call it lot 3B. (Shirley and David have 3A, more concentrated.) This stuff is at 10<sup>4</sup> $\mu$ l using their activated salmon sperm (or is it calf thymus) DNA assay. Titrate for PCR amplification.

A, F: 1 $\mu$ l per 100 $\mu$ l rxn.  
B, G: 1/2 "  
C, H: 1/4 "  
D, I: 1/8 "  
E, J: 1/16 "

A-E: Molt4  
F-J: GM2064

Molt4 and GM2064 @ 100 $\mu$ g/ $\mu$ l  
PC03 and PC04 @ 10 $\mu$ M, dNTP @ 10mM each

35 $\mu$ l Molt4 or GM2064  
35 $\mu$ l 10x Tag salts  
35 $\mu$ l PC03  
35 $\mu$ l PC04  
35 $\mu$ l DMSO  
52.5 $\mu$ l dNTP  
122.5 $\mu$ l H<sub>2</sub>O  
350 $\mu$ l  $\rightarrow$  5', 95°

Cool to RT and divide into one 100 $\mu$ l <sup>sample</sup> ~~volume~~ and four 50 $\mu$ l samples (50 $\mu$ l left over). Add 1 $\mu$ l Tag polymerase (lot 3B, 10<sup>4</sup> $\mu$ l) to the 100 $\mu$ l volume and mix. Prepare four 50 $\mu$ l serial dilutions in four remaining samples. <sub>two-fold</sub>

Final concentrations of enzyme per 100 $\mu$ l reaction volume: 1 $\mu$ l, 1/2, 1/4, 1/8, 1/16.

Overlay with mineral oil

To Page No. 82

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

*[Signature]*

*R. Sankar*

From Page No 81


Save remaining 50  $\mu$ l ( $\frac{1}{16}$  dil'n) and store @ 4° (just in case enzyme works at  $< \frac{1}{16}$   $\mu$ l per 100  $\mu$ l).

Subject to 24 cycles: 2 min ramp,  $35^{\circ}$  to  $95^{\circ}$   
2 min ramp,  $95^{\circ}$  to  $35^{\circ}$

After last cycle, incubate additional 5 min at  $65^{\circ}$  to complete final (25<sup>th</sup>) extension.

Extract oil with  $\text{CHCl}_3$ .

Load 5 $\mu$ l each A-D and F-I on 4% NuSieve/0.5% agarose/  
1xTBE.



ШОПКАМ ГАМКАМ  
НА УРАКАМ

Got PCR product in Molt4  
at all four dilutions,  
even  $1/8 \mu\text{l}$ ! As expected,  
nothing in GM2064.

Background is very low, virtually nonexistent for  $1/4$  and  $1/8 \mu\text{l}$  samples. Maybe combination of enzyme and "fast ramp" protocol is responsible.

Need to check 1/16 ul samples; might be band there, too.

To Page No. 83

Witnessed & Understood by me,

Date \_\_\_\_\_

**Invented by**

Date \_\_\_\_\_

**Recorded by**

Compare titrations of lot 2A and lot 3B of Tag polymerase using "Pro/Pette" protocol.

A I: 1  $\mu$ l  
B J:  $\frac{1}{2}$   
C K:  $\frac{1}{4}$   
D L:  $\frac{1}{8}$

E M:  $\frac{1}{16}$   $\mu$ l  
F N:  $\frac{1}{32}$   
G O:  $\frac{1}{64}$   
H P:  $\frac{1}{128}$

A-H: lot 2A  
I-P: lot 3B

50  $\mu$ l Melt4  
50  $\mu$ l 10x salts  
50  $\mu$ l PC03  
50  $\mu$ l PC04  
50  $\mu$ l DMSO  
75  $\mu$ l dNTP  
175  $\mu$ l H<sub>2</sub>O  
500  $\mu$ l  $\rightarrow$  10', 95°

Prepare eight 50  $\mu$ l two-fold serial dilutions with 1  $\mu$ l lot 2A or lot 3B as described on page 88.

Subject to 24 cycles in ProPette: 2 1/2' min ramp, 37 to 95°  
3' min ramp, 95 to 37°

After last cycle incubate 5' @ 60° to complete 25<sup>th</sup> cycle extension.

Load 5  $\mu$ l each onto 4% NuSieve/0.5% agarose/1x TBE.

To Page No. 92

Read & Understood by me,

Date

Invented by

Date


Recorded by

*Elaine*

*R. Saito*

No. 90

A B C D E F G H I J K L M N O P



amplification

Don't see cut-off in lot 2A anymore. Can see a band as far down as  $1/64 \mu\text{l}$  (lane G). Also, best S/N ratio is at  $1/16 \mu\text{l}$  instead of  $1/8 \mu\text{l}$ .

lot 3B seems to peak at  $1/2 \mu\text{l}$  ~~but~~ although this should be rechecked. ~~(Especially recheck photo on page 1049 shows peak at 1/2 and 1/4)~~

There seems to be more background in these samples than in those done "fast ramp". ~~Compare~~ Need to compare on same gel to be sure.

This gel show effect of polymerase on S/N ratio quite nicely.

To Page No. X

Read &amp; Understood by me;

Date

Invented by

Date

Recorded by

R. Saiti

From Page No. X

Lot 3 Tag polymerase seems to be losing activity. Initially, optimal concentration (for PC03/04) was at 2.5 u per 100  $\mu$ l reaction. Over the last 2-3 weeks activity has gradually ~~disappear~~ dissipated. Most recent attempts have failed (not recorded). Only a very weak PCR ~~products~~ product was seen with 5u. Russ, Dory, and Steve have had similar experiences.

Unlike lot 2A, the storage buffer for lot 3 does not contain the non-ionic detergents Tween 20 or NP40.\* David Gelfand's experience with this polymerase indicates that it is a sticky enzyme and he routinely uses both detergents during purification to improve yield and during assay to stimulate activity.

May be that in the absence of detergents and at -20° the enzyme is aggregating. Addition of "soap" to either the storage buffer or the PCR reaction may restore activity. Will try the latter first.

	$\mu$ /100 $\mu$ l
A H:	$\frac{1}{2}$
B I:	$\frac{1}{4}$
C J:	$\frac{1}{8}$
D K:	$\frac{1}{16}$
E L:	$\frac{1}{32}$
F M:	$\frac{1}{64}$
G N:	$\frac{1}{128}$

A-G: (-) detergent  
H-N: (+) detergent (0.05% each)

\* Another difference is that lot 3 contains 200  $\mu$ g/ml gelatin. Lot 2 doesn't.

To Page No. 102

Witnessed & Understood by me,

*John F. ...*

Date

Invented by


*R. Santi*

Date

Recorded by

Load 5 $\mu$ l each onto 4% NuSieve/0.5% agarose/1xTBE  
→ 100v, 90'

A B C D E F G H I J K L M N



Detergents definitely have some effect. Without them can only see a band in ~~the~~ 5 $\mu$  sample (A). But with them, can see bands in 5 $\mu$  (H), 2.5 (I), and very faintly in 1.25 (K).

Although Tween and NP40 helped, activity still is not as good as it was originally (p. 81).

Maybe that adding more detergents or adding instead to the storage buffer will ~~that~~ be better.

Molt4 @ 100  $\mu\text{g}/\text{ml}$ , PC03 and PC04 @ ~~10  $\mu\text{M}$~~  10  $\mu\text{M}$   
 dNTP @ 40 mM (p.99), Tween/NP40 @ 0.5% each

A-G:	50 $\mu\text{l}$	Molt4	H-N:	50 $\mu\text{l}$	Molt4
	50 $\mu\text{l}$	10x salts		50 $\mu\text{l}$	10x
	50 $\mu\text{l}$	PC03		50 $\mu\text{l}$	PC03
	50 $\mu\text{l}$	PC04		50 $\mu\text{l}$	PC04
	50 $\mu\text{l}$	DMSO		50 $\mu\text{l}$	DMSO
	75 $\mu\text{l}$	dNTP		75 $\mu\text{l}$	dNTP
	175 $\mu\text{l}$	H <sub>2</sub> O		50 $\mu\text{l}$	Tween/NP40
	500 $\mu\text{l}$	$\rightarrow 10', 95^\circ$		125 $\mu\text{l}$	H <sub>2</sub> O
				500 $\mu\text{l}$	$\rightarrow 10', 95^\circ$

Dispense each 500  $\mu\text{l}$  mix into one 100  $\mu\text{l}$  sample and six 50  $\mu\text{l}$  samples. Add  $1/2 \mu\text{l}$  lot 3 to the 100  $\mu\text{l}$  sample and ~~dilute~~ serially dilute, preparing seven 50  $\mu\text{l}$  two-fold serial dilutions:  $1/2 \mu\text{l}$  to  $1/128 \mu\text{l}$ .

Overlay with mineral oil and subject to 24 cycles on Sinsky's ProPette: 3' ramp  $37^\circ$  to  $95^\circ$  (hot water set at  $102^\circ$ )  
 3' ramp  $95^\circ$  to  $37^\circ$

After last cycle, incubate additional 10' at  $56^\circ$ .

Extract oil with  $\text{CHCl}_3$ . (Samples H-N ~~became~~ became cloudy. Probably interaction of detergents with chloroform.)



From Page No. X

Follow up on expt. described page 101 by ~~determining~~ seeing if adding detergents directly to enzyme stock is a better way to go. Gelfand suggests 0.5% each is a good starting point.

19  $\mu$ l Tag polymerase (lot 3,  $10^4$   $\mu$ l)  
 1  $\mu$ l 10% Tween 20 / 10% NP-40  
 20  $\mu$ l Tag pol,  $9.5^4$   $\mu$ l, 0.5% each detergent  $\rightarrow$  incubate @ RT for ~10', mixing thoroughly. Store @  $4^\circ$

A E:	1	} $\mu$ l Tag	A-D:	Tag w/o detergent
B F:	$\frac{1}{2}$		E-H:	Tag w/ detergent
C G:	$\frac{1}{4}$			
D H:	$\frac{1}{8}$			

reagents as desc. page 101

30  $\mu$ l Molt4  
 30  $\mu$ l 10x salts  
 30  $\mu$ l PC03  
 30  $\mu$ l PC04  
 30  $\mu$ l DMSO  
 45  $\mu$ l dNTP  
 105  $\mu$ l H<sub>2</sub>O  
 300  $\mu$ l  $\rightarrow$  10',  $95^\circ$

Prepare two 300  $\mu$ l mixes. Divide each into one 100  $\mu$ l sample and three 50  $\mu$ l samples. Add 1.0  $\mu$ l enzyme, with or without detergent, to the 100  $\mu$ l sample and serially dilute 50  $\mu$ l into the three 50  $\mu$ l samples.

Amplify and workup as described p102 except use our Pro/Pette with this program:  $2\frac{1}{2}$  ramp,  $35^\circ$  to  $95^\circ$   
 3' ramp,  $95^\circ$  to  $35^\circ$

Ward Smith

To Page No. 105

Witnessed &amp; Understood by me,

Date

Invented by

Date


Recorded by

Tag PCR: (cont'd)

Page No. 104

Load 5 $\mu$ l each sample onto 4% NuSieve/0.5% agarose/1x TBE.  
 → 100V

A B C D E F G H



This is it! Activity of enzyme with detergent is as good as (maybe even better) than original titration (see p. 81).

Based on this expt. best conc. of enzyme is either 5.0 u (F) or 2.5 u (G). Former may have a teeny bit more PCR product, but latter has less background. (Either is fuckin' good.)

Looks as if activity in (-) detergent enzyme has gotten even worse. Can barely see the ~~4th~~ 1/2 $\mu$ l sample (B).

Should add detergents to the remaining enzyme stock.

To Page No. X

Witnessed &amp; Understood by me,

Date

Invented by

Date

Recorded by

*Jul Fabian*

*R. Saini*